



DOCTORS



ECO Gel

**SAF DOCTORS
ECO GEL**



Why Choosing SAF Doctors Eco Gel?


SAF Doctors ECO gel is suitable for both ultrasound and laser hair removal (IPL) procedures :

- Skin-friendly harmless to sound heads
- Microbubble-free,
- Can be spread very thinly and
- Will not dry out prematurely,
- Water soluble
- Will not stain clothing
- Does not contain harmful chemicals,
- Free from formaldehyde and alcohol.



How to use:

Add sufficient amount on the area and spread the gel for usage.



Types of SAE Doctors ECO Gel:

SAF Doctors Ultrasound gel 250 ml bottle is easy to dispense gel from and easy to handle.

Its dispenser designed to stop dripping.

SAF Doctors Ultrasound Gel is the mostly used gel for examination ultrasound and IPL.

Can be both color free and with color.



SAF Doctors ultrasound 500 ml bottle is easy to dispense gel from and easy to handle.

Its dispenser designed to stop dripping.

The 500 ml useful for operations where more gel is required, like sonographic procedure, IPL and more.

Can be both color free and with color.



SAF Doctors Gallon Ultrasound Gel is the most cost-effective way to purchase gel.

Each Gallon comes with an empty 250 ml dispenser bottle for easy dispensing.



SAF Doctors 1000 ml Nylon Pack:

is the most economical way to purchase Eco gel, is designed in way where every drop in the pack to be used



Products Codes:

- 1000008: ECO GEL GAL Transparent ≈ 5 L
- 1000009: ECO GEL GAL. blue ≈ 5 L
- 1000010: ECO GEL Transparent bott. 500ml-10pcs/box
- 1000011: ECO GEL blue bott.500-10pcs/box
- 1000012: ECO GEL Transparent bott. 250ml-12pcs/box
- 1000013: ECO GEL blue bott. 250-12pcs/box
- 1000002:ECO GEL 1000 ML NYLON PACK